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
Detection Results				QC Information			
1. List of Bacteria Detected				Item	QC Parameter	Detection Results	QC Standard
Type	Name	Latin Name	Seq	Identification of pathogen	Positive control (PC)	Qualified	Qualified
G-	<i>Prevotella melaninogenica</i>	<i>Prevotella melaninogenica</i>		Negative control (NC)	Qualified	Qualified	
Type G (Gram-positive bacteria) / G (Gram-negative bacteria)				Anticopy number	Detection rate of internal reference	100%	100%
				Overall evaluation of sequencing QC	Pass	Qualified	
2. List of Fungi Detected				1. Negative control (NC): Analyze the contamination of the whole detection process; 2. Artificial sequence is selected as PC, and the detection of PC can reflect whether the test is successful.			
7. List of Other Pathogen Detected				3. List of Drug Resistance Genes Detected			
Type	Name	Latin Name	Number of Sequences Detected	Relative	Gene Name	Resistant Drug Category	Gene Information
MC	<i>Chlamydia psittaci</i>	<i>Chlamydia psittaci</i>	25		WZO	tetracycline antibiotics	WZO is a ribosomal protection protein that confers resistance to tetracycline antibiotics. Its genes are associated with genetic transposons and are found in both Gram-positive and Gram-negative bacteria. Its resistance mechanism is mainly to protect the target of antibiotic action from the influence of antibiotic combination.
8. List of Suspected Background Microorganisms							
Type	Name	Latin Name	Number of Sequences Detected	Relative			
G+	<i>Neisseria subflava</i>	<i>Neisseria subflava</i>	400	10.00%			
G+	<i>Streptococcus sanguinis</i>	<i>Streptococcus sanguinis</i>	90	2.25%			
G-	<i>Haemophilus parvulus</i>	<i>Haemophilus parvulus</i>	31	0.77%			
G-	<i>Corynebacterium granulosum</i>	<i>Corynebacterium granulosum</i>	15	0.37%			

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1. Paul J. Introduction to Infectious Diseases[M]//Disease Causing Microbes. Cham: Springer International Publishing, 2024: 1-63.
2. Huang Z, Wu Q, Fang X, et al. Comparison of culture and broad-range polymerase chain reaction methods for diagnosing periprosthetic joint infection: analysis of joint fluid, periprosthetic tissue, and sonicated fluid[J]. International orthopaedics, 2018, 42: 2035-2040.
3. McCombie W R, McPherson J D, Mardis E R. Next-generation sequencing technologies[J]. Cold Spring Harbor perspectives in medicine, 2019, 9(11): a036798.
4. Chiu C Y, Miller S A. Clinical metagenomics[J]. Nature Reviews Genetics, 2019, 20(6): 341-355.
5. Li S, Tong J, Li H, et al. L. pneumophila Infection Diagnosed by tNGS in a Lady with Lymphadenopathy[J]. Infection and Drug Resistance, 2023: 4435-4442.
6. Chao L, Li J, Zhang Y, et al. Application of next generation sequencing-based rapid detection platform for microbiological diagnosis and drug resistance prediction in acute lower respiratory infection[J]. Annals of Translational Medicine, 2020, 8(24).
7. Gao D, Hu Y, Jiang X, et al. Applying the pathogen-targeted next-generation sequencing method to pathogen identification in cerebrospinal fluid[J]. Annals of Translational Medicine, 2021, 9(22).
8. Huang C, Chen H, Ding Y, et al. A microbial world: could metagenomic next-generation sequencing be involved in acute respiratory failure?[J]. Frontiers in cellular and infection microbiology, 2021, 11: 738074.
9. Li B, Xu L, Guo Q, et al. GenSeizer: a multiplex PCR-based targeted gene sequencing platform for rapid and accurate identification of major mycobacterium species[J]. Journal of Clinical Microbiology, 2021, 59(2): 10.1128/jcm. 00584-20.

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